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
Prof. Emanuele Montomoli Tel +39/0577/234134 e-mail: montomoli@unisi.it

Final report:

MICROBIOLOGICAL CONTROLS EFFECTUATED ON THE AIR IN A ROOM IN WHOM IS COLLOCATED AN OZONISED HYDRO MASSAGE SYSTEM PRODUCED BY OZONOMATIC  $^{\tiny \odot}$ 

	Firma	Data
Responsabile del Laboratorio	Elmontant	20-Feb-2004
Esecutore esperimento	Chory Retile	20-Feb-2004

Oranformed of Repositions, Millians
Sportmormers and Paulician
Sportmormers and Paulician
Universitä dengli Studiel of Sportmormers
Link Epidemiologia Milliane langua
America Paul Zanasa ali dan barana
America Paul Zanasa ali dan barana

Microbiological controls effectuated on the air/ water of aonized hydro massage system produced by Ozonomatic

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Prof. Emanuele Montomoli Tel +39/0577/234134 e-mail: montomoli@unisi.it

### Final report:

# MICROBIOLOGICAL CONTROLS EFFECTUATED ON THE WATER OF AN OZONISED HYDRO MASSAGE SYSTEM PRODUCED BY OZONOMATIC®

	Firma	Data
Responsabile del Laboratorio	Elmontant	20-Feb-2004
Esecutore esperimento	Chord Rettle	20-Feb-2004

Microbiological controls effectuated on the water of an ozonized hydro massage system produced by Ozonomatic

### Role of hydro-ozono in the treatment of diabetic foot.

M. Martinelli, S. Rotunno, A. Armiento, A. Bianchi, O Guarino, L. Giubilei, F. Lasaracina, C. Valente, L. Xiao, M. Cassol Divisione di Medicina Interna e Geriatria- Ospedale San Pietro, Fatebenefratelli-Rome

### Background:

The World Health Organization has estimated 360 million of diabetics in 2012 compared with 120 million estimated in 1996; Other OMS estimates show that about 15% of diabetics, are subject to the foot ulcer which requires major medical treatment (54 million patients estimated to 2012). The biggest problem related to an ulcer in the diabetic foot is the risk of partial or total amputation of the foot. Diabetics provide more than 50% of all amputations from the ankle and up.

The vascular disease of the diabetic foot is essentially a microangiopathy, although large vessels have multiple calcified stenosis all along their course. These patients belong to a subgroup of diabetics with retinopathy-nephropathy-peripheral microangiopathy, clinical aspects of the "diabetic capillaritis" that emerges by the appearance of microalbuminuria and proteinuria in the urine.

The effect of microangiopathy, reduced capacity to resist caused by to repeated microtrauma capillaries caused by iperpression; reduced capacity to repair lesions with neovascularization. The result is the increased facilities onset of trophic lesions, their tendency to slow healing and gangrene.

### Am of the study:

We have carried out a study to verify the antibacterial effect of ozone and if it could be applied in disinfection of diabetic foot ulcers. Ozone when added to water reacts by forming some types of reactive oxygen species (ROS) which have disinfectant activity. (Fig.1)

### Methods

Valutation of the effect on the bacterial charge present in the water following up the imput of ozone in a foot bath; his was performed by using an ozonized whirlpool (Ozonomatic System-OS) (Fig.2). It works by using oxygen in the air and converting into very safe levels of ozono by special lamps (Corona effect).

A mixed cultivation ground (E.Coli, S.Aureus, P. Aeruginosa, S. Faecalis and L.Pneumofila) has been added to 201 of water (T 37°C).

After shaking, two samples has been withdrawal of 10 ml and 5 ml of wather ("T0") and have been incubated in slabs with cultivation ground appropriated for every five bacterial part. At the and of the cycle of emission of ozone in the water (20"), time "T1", have been repeated with the same methods as T0 (Tab. 1).

### Discussion:

From the evaluation of the results has come out that the ozonized whirlpool system by Ozonomaticsystem offer an effective diminution of the microbial charge after the 20 minutes stabilised for the ozone effect. The strength of disinfection can be increased in two ways:

- 1) increasing to 30 minutes the bath-ozonized time
- 2) by replacing or adding to the drinking water the sterile saline.

The ROS forming whith the eposure of water to ozone reacts with NaCl (contrary to the findings in the literature) producing:

Sodium hypochlorite - NaOCI

Hypochlorous acid - HOCI,


Their radical - OCI (highly reactive molecule with high disinfection power).

The immersion of the foot in the water disinfectant, at 37  $^{\circ}$  C for 30 minutes, results in:

- increase in the period of contact between bacteria and disinfectant that is also continuously regenerated
- the wound is softened and made more permeable to disinfectants that also reach the deeper layers of tissue as result of the vasodilation due to the temperature of the bath
- the possibility that the disinfectant is able to enter even in the surface layers of exposed bone. (Tab.2; Fig.3)

### Conclusion:


The ozonized whirlpool system by Ozonomaticsystem has been shown to have excellent ability to reduce the bacterial charge and disinfect the diabetic foot ulcers.

	T0 UFC/10ml Average	T0 UFC/5ml Average	T1 UFC/10ml Average	T1 UFC/5ml Average
Bacteria total at 36 °C	1392	812	380	263
Escherichia coli	812	464	580	348
Pseudomonas aeruginosa	1856	1508	738	668
Strafilococcus aureus	2900	2320	35	25
Streptococcus faecalis	1560	840	464	348

### Table1:

T0; water sample before Ozone treatment T1; water sample after Ozone treatment

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### Table1:

T0; water sample before Ozone treatment T1; water sample after Ozone treatment

Figure1: hydroozone therapy can acts

Figure2: Ozonomatic System


Table2:REDUCTION (%) of the bacteria

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EC Quality Assurance System Certificate

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Via F. P. Bonifacio, 52/54 00156 Roma (RM)

Italia |

Exertises simple to applicable to the This completion with the applicable equipment of

Direttiva 93/42/CEE, Allegato V, attuata in Italia con Dlgs. 46 del 1997/02/24 93/42/EEC Directive, Annex V, transposed in Buly by Digs. 46 of 1997/02/24

per la request springe di Disposites Medici i de diviniery Michael Zirecus.

Idroterapia Ozonizzata / Ozonized hydrotherapy

Month Comment of Contract Cont

Organismo Notificato n. 0476

R provide Certificate a suggests at reports del regissements CERMET of a wallet take per la typologie di dispositie supra abrello dei suggests a conseptiones. L'adrepte secretar parte receptaire del presente Certificate.

The Cartholic Endoctor CLEMIT registrates and it is sold only be the observe mentioned Minday Devices that are subject to somey. The soldness areas is an imprompt part of

Date of professional hours duty Data of time modifical Last strategistation Data wardening Controller 20579030

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DGFDMIIII.S.Le.1/ /P.

HAVING REGARD to the 93/42/EEC Directive concerning medical devices;

HAVING REGARD to the Legislative Decree (D.L.vo) no. 46/97 (and its following amendments) reporting the accomplishment of directive 93/42/EEC:

HAVING REGARD to the request dated 15 January 2009, reference nr. 4244-A-04/02/2009, submitted by the Company OZONOMATIC SYSTEM S.F.I., located in 00156 Roma (RM) Italy, Via F.P. Bonifacio 52/54.Vat. N. IT06342111009:

WHEREAS this Company paid the fees required by Ministerial Decree (D.M.) May 24, 2004; MAVING REGARD to the official deeds:

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SIATTESTA

That the Company OZONOMATIC SYSTEM S.r.I., located in 00156 Roma (RM) Baly , Via F.P. Bonfacio 52/54, Vat. N. IT06342111009: has manufactured and marked CE as medical device. according to the Directive 93/42/EEC, the following: product:

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The above mentioned product, according to the art. 4 of the above mentioned Directive, can freely circulate and can be placed on the market in Italy and in all the European Union.

This certificate is issued on the interested party's request in according to the law and for export to Countries outside European Union.

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

RELATORE: Dott. Diego Ceci

A.A. 2015/2016

### Role of hydro-ozono in the treatment of diabetic foot.

M. Martinelli, S. Rotunno, A. Armiento, A. Bianchi, O Guarino, L. Giubilei, F. Lasaracina, C. Valente, L. Xiao, M. Cassol Divisione di Medicina Interna e Geriatria- Ospedale San Pietro, Fatebenefratelli-Rome

### Background:

The World Health Organization has estimated 360 million of diabetics in 2012 compared with 120 million estimated in 1996; Other OMS estimates show that about 15% of diabetics, are subject to the foot ulcer which requires major medical treatment (54 million patients estimated to 2012). The biggest problem related to an ulcer in the diabetic foot is the risk of partial or total amputation of the foot.

Diabetics provide more than 50% of all amputations from the ankle and up.

The vascular disease of the diabetic foot is essentially a microangiopathy, although large vessels have multiple calcified stenosis all along their course. These patients belong to a subgroup of diabetics with retinopathy-nephropathy-peripheral microangiopathy, clinical aspects of the "diabetic capillaritis" that emerges by the appearance of microalbuminuria and proteinuria in the urine.

The effect of microangiopathy, reduced capacity to resist caused by to repeated microtrauma capillaries caused by iperpression; reduced capacity to repair lesions with neovascularization. The result is the increased facilities onset of trophic lesions, their tendency to slow healing and gangrene.

### Am of the study:

We have carried out a study to verify the antibacterial effect of ozone and if it could be applied in disinfection of diabetic foot ulcers. Ozone when added to water reacts by forming some types of reactive oxygen species (ROS) which have disinfectant activity. (Fig.1)

### Methods:

Valutation of the effect on the bacterial charge present in the water following up the imput of ozone in a foot bath; his was performed by using an ozonized whirlpool (Ozonomatic System-OS) (Fig.2). It works by using oxygen in the air and converting into very safe levels of ozono by special lamps (Corona effect).

A mixed cultivation ground (E.Coli, S.Aureus, P. Aeruginosa, S. Faecalis and L.Pneumofila) has been added to 20 I of water (T 37°C).

After shaking, two samples has been withdrawal of 10 ml and 5 ml of wather ("T0") and have been incubated in slabs with cultivation ground appropriated for every five bacterial part. At the and of the cycle of emission of ozone in the water (20"), time "T1", have been repeated with the same methods as T0 (Tab. 1).

### Discussion:

From the evaluation of the results has come out that the ozonized whirlpool system by Ozonomaticsystem offer an effective diminution of the microbial charge after the 20 minutes stabilised for the ozone effect. The strength of disinfection can be increased in two ways:

- 1) increasing to 30 minutes the bath-ozonized time
- 2) by replacing or adding to the drinking water the sterile saline.
- The ROS forming whith the eposure of water to ozone reacts with NaCl (contrary to the findings in the literature) producing:

Sodium hypochlorite - NaOCI

Hypochlorous acid - HOCI,

Their radical - OCI (highly reactive molecule with high disinfection

The immersion of the foot in the water disinfectant, at 37 ° C for 30 minutes, results in:

- 1) increase in the period of contact between bacteria and disinfectant that is also continuously regenerated
- 2) the wound is softened and made more permeable to disinfectants that also reach the deeper layers of tissue as result of the vasodilation due to the temperature of the bath
- 3) the possibility that the disinfectant is able to enter even in the surface layers of exposed bone.(Tab.2; Fig.3)

The ozonized whirlpool system by Ozonomaticsystem has been shown to have excellent ability to reduce the bacterial charge and disinfect the diabetic foot ulcers.

NEUROPATHY	MICRO ANGSOPATHY WASCULOPATHY
REDUCED PAIN SENSITIVITY, REDUCED PAIN ALARM  CUTANEOUS HYPER PRESSION	REDUCED CAPILLARY RESISTANCE TO MICRO TRAINA IMPAIRED ARILITY
OZONE GANG	ING
АНРИЛ	INFECTION

T0 UFC/10ml Average	T0 UFC/5ml Average	T1 UFC/10ml Average	T1 UFC/5ml Average
1392	812	380	263
812	464	580	348
1856	1508	738	668
2900	2320	35	25
1560	840	464	348
	UFC/10ml Average 1392 812 1856 2900	UFC/10ml UFC/5ml Average Average 1392 812 812 464 1856 1508 2900 2320	UFC/10ml Average         UFC/5ml Average         UFC/10ml Average           1392         812         380           812         464         580           1856         1508         738           2900         2320         35

T0; water sample before Ozone treatmen T1: water sample after Ozone trea



Table2:REDUCTION (%) of the bacteria

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Biomed Res Int. 2017;2017:5842723.

doi: 10.1155/2017/5842723. Epub 2017 Mar 2.

Thermography Improves Clinical Assessment in Patients with Systemic Sclerosis Treated with Ozone Therapy.

Nowicka D1.

Department of Dermatology, Venereology and Allergology, Wrocław Medical University, Wrocław, Poland.

### **Abstract**

Objective. Treatment of scleroderma is challenging and limited. The aim of our study was to evaluate the usefulness of thermography in assessment of the clinical condition (joints movability and skin thickness) in clinically advanced patients with systemic sclerosis before and after ozone therapy.

Method. The study included 42 patients aged 32 to 73 years with advanced systemic sclerosis hospitalized in the university clinic between 2003 and 2006. Thermography and clinical examinations were conducted at baseline and after two series of bath in water with ozone.

Results. The comparison of results showed significant increase in skin temperature by  $2.5^{\circ}$ C, significant increase in interphalangeal joints movability by 18 degrees, and significant decrease in skin score by 14.7 points. The skin temperature was correlated with skin score (r = -0.59) and joints movability (r = +0.8).


**Conclusions.** Ozone therapy shows positive effect on clinical parameters and skin temperature as measured with thermography. The study indicated possibility of introducing ozonotherapy as an independent therapy in cases with low level of progression or during remission periods and as additional treatment in patients with advanced disease requiring immunosuppressive treatment. Thermography is useful in assessment of skin condition showing strong correlation between skin temperature and clinical parameters.

PMID: 28349063 PMCID: PMC5352863 DOI: 10.1155/2017/5842723

[Indexed for MEDLINE] Free PMC Article

Thermographic pictures with corresponding histogram of color scale.

- (a) A healthy hand: the noticeably lighter places in the area of distal finger parts and the area of the ball point to higher skin temperature in those places.
- (b) A hand of a person suffering from systemic sclerosis before ozone therapy: significant lowering of the temperature of phalanx distal parts (especially of the 4th finger) is visible.
- (c) A hand of a person suffering from scleroderma after ozonotherapy treatments: the increase in the temperature in the entire examined area manifesting itself as lighter colors on the thermogram is visible.

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